

# hospicall P7

• Nurse Call Light Systems

A

- Personal Protection Systems
- Emergency Call Systems

... just press the button – and help is on the way!

**Dedicated, sustained, dynamic.** This is our mission statement, this is hospicall GmbH. We have our roots in Wiehl, in the picturesque Oberberg district east of Cologne. This is where we develop, produce, and live.

For more than 15 years we have been meeting the challenges of the market – producing call systems that set standards both economically and innovatively. Our staff's experience and know-how are key factors to our success.

Far beyond our core capabilities as developers and producers, we see ourselves as partners passionately working every single day to provide optimal and flexible solutions that are ideally tailored to a customer's specific applications and requirements. Our objective is to facilitate product installation, and to provide reliable safety for the end user. Just press the button – and help is on the way! To cut a long story short – stateof-the-art technology and reliability both in the clinical environment and for personal protection are our motivation and incentive.

This in mind, we are continually refining our products step by step to present innovative and enhanced solutions at any given time – now and also in the future. This is what makes us authentic, this is what hospicall GmbH stands for.

# Dedicated. Sustained. Dynamic.



System Description

P7 System Overview (with

P7 System Overview (witho

P7 System Control Unit, Mo

FAP Interface, DECT Interfa DCF77 ISO Interface, 19" Bi

SBUS Splitter, RS232 Isolat

Mains Adapter 24V 5A, Main

Power Rectifier with UPS M Power Rectifier with UPS M

Ward Controller P7, Standa

Galvanic Ethernet Isolator,

Duty Room Monitor

Room Terminal Duty Room Table Top Stand for Duty Ro

Room Terminal with Displa

Room Modules, Room Mod

Call Units, Display Module

Presence / Cancel Push-Bu

Call Units, Pull Cord Switch

Call Units, Diagnostic Unit,

Multiple Push-Button – Ton Wall Bracket for Multiple Pu with Auxiliary Socket and D

Call Unit with Auxiliary Plu Plug Contact and Tone, Cal Trilux Bedside Light, Call U

Patient Handset, Wall Brac Wall Bracket for Pear / Mult

LED Corridor Indicator lam Nameplate with LED Corric a Corridor Indicator lamp, I

Corridor Display, Single-Sid Double-Sided for Fixation t

Interface Door Intercom, Co Output Module, I/O-2AT Inp

I/O-8 Input / Output Modul

Radio Sets

Tone Detector

Sealing Flange for fm Inser

Push-Button Program





www.hospicall.com



|   | 4 - 5                      |
|---|----------------------------|
| Speech Function)  | 6 - 7                      |
| out Speech Function)  | 8                          |
| obile Call Display, SIP-Gateway   | 10                         |
| ice, Graphic Message Display,<br>racket, Isolator   | 11                         |
| or 5kV  | 12                         |
| ns Adapter 24V 10A  | 14                         |
| Iodule 24V, 5A, incl. Rechargeable Battery<br>Iodule 24V, 10A, incl. Rechargeable Battery   | 15                         |
| alone Ward Controller, Serial Interface   | 16                         |
| Ethernet Switch 8 Port, Ethernet Switch 5 Port  | 17                         |
|   | 18                         |
| n with Display, with/without Speech Function,<br>oom Terminal   | 20-21                      |
| ay, with/without Speech Function  | 22-23                      |
| lule with Blank Cover   | 24-25                      |
|   | 26-27                      |
| ttons, Reader Module, Transponder for Room Electronics  | 28-29                      |
| n, Pneumatic Call Button, Replacement Pull Cord   | 30-31                      |
| Loudspeaker Unit, Audio Isolator  | 32-33                      |
| ne, incl. Self-releasing Protective Adapter,<br>ush-Button – Tone, Call / Cancel Push-Button<br>Diagnostic Socket   | 34                         |
| g Contact, Diagnostic Socket and Tone, Call Unit with Auxili<br>Il Unit with Auxiliary Plug Contact, Diagnostic Socket and To<br>Init with Auxiliary Plug Contact and Tone for Trilux Bedside | ary<br>one for<br>Light 35 |
| ket, Connection Extension, Pear / Multiple Push-Button,<br>tiple Push-Button, Trapeze Bar Clip  | 36-37                      |
| p, Doorplate with LED Corridor Indicator lamp,<br>dor Indicator lamp, Electronics for connecting<br>Nameplate Service Set   | 38-39                      |
| ded for Fixation to Wall, Corridor Display,<br>to Ceiling   | 40                         |
| urrent Impulse Switch, I/O Input /<br>put / Output Module with Cancel Push-Button   | 42-43                      |
| le  | 44                         |
|   | 46-47                      |
|   | 48                         |
| ts, BUS Wall Socket   | 49                         |
|   | 50-51                      |

# **System Description**

#### General

The call system hospicall P7 is a modern call system which conforms to DIN VDE 0834 1/2:2000-04 and takes account of, in particular, the increasing economic requirements of modern call systems. This applies not only to the purchase and installation but also to the operation in the day-to-day care environment of hospitals and care homes. The system can be used in the Application Areas 1 and 2 according to DIN / VDE.

#### System Structure

hospicall P7 is a de-centrally structured, scalable system, which, through its modularity and flexibility, is just as suitable for the smallest systems with two or three treatment rooms as for an institution with several different residential areas or wards and function areas.

Hereby, the System Control Unit is the master control unit for the entire system and, in addition to the required monitoring functions, takes on responsibility for communications spanning the zones between the connected zone controllers and third-party systems. The ward controller serves as the control unit within a ward and regulates the data traffic between the components in the room and corridor.

An important role of the control is the optical and acoustic call forwarding and/or call handling as well as the forwarding of this information to the higher ranking modules. The room terminals, room modules and information displays are connected to the ward controller via the ward bus. In particular, the room terminals and room modules act as control units for a room. As well as having large controls for the care staff, the room terminal has a large, easily readable colour touch-display for indicating and processing calls and the necessary components for speech communication, such as loudspeaker and microphone.

Depending on the how the room is equipped, the other release elements of the room as well as the corridor indicator lamps (where not integrated into the room module) are connected to the room terminal or room module. Different call units can be connected for the various call release requirements. In the various models, the call units are supplied with large, easily operated call and cancel pushbuttons, as well as presence push-buttons with membrane cover for easy cleaning and disinfecting. Beneath all membrane-covered operating surfaces there are micro switches, which provide a good tactile feel.

Jack sockets in the call units enable the connection of pear and multiple pushbuttons and patient handsets. They can also be used for the connection of medical equipment, which can send alarms via the call system.

In the corridor the care personnel are informed via the large wall or ceiling mounted displays; duty room terminals for wall mounting or desktop table mounting are available in the duty room, by means of which organisational measures such as zone linking and general announcements can be handled, as well as call processing.

#### System Integration / Interfaces

Of course the system hospicall P7 has a variety of interfaces available for messaging systems (DECT systems, paging systems, GSM mobile telephone) or old equipment. The call logging and evaluation can be made available in each duty room via a standard data network (CAT5 wiring according to DIN EN 50173). The operating devices in the rooms with their mechanical characteristics are adapted to the current requirements of the switch industry, taking into account the architectural aspects, which are not only indispensable in pleasantly laid out senior residences. Also, the integration of these units in lighting rails or medical supply units is guaranteed.

#### Configuration

The system configures itself on commissioning. All room and ward bus clients are automatically recognised and are provided with function typical parameters (WC units, bed units).

The configuration and setting up of features, e.g. room-by-room reduction in volume of call forwarding, takes place from a central position, independently of the operating system, with a browser operator interface.











#### Supply Network

For the wiring of the call system hospicall P7, the telecommunication cable type IY(St) Y is sufficient. Special system cables are not required. This means that, above all, the use of existing cables when renovating is a great economic advantage. Apart from the wires for the current supply (24V, GND), the data lines and speech lines each require a twisted pair of wires. The length of the wire can thereby be up to 1000 m.

#### Safety

Since with a call system one is dealing with safety equipment, the system hospicall P7 is fitted with different safety functions. On the one hand, short-circuit-proof inputs protect the electronics if wiring errors occur in the installation phase; on the other hand, the central hierarchical system structure guarantees the basic functions of the release and indication of calls even if individual modules fail.

As the smallest unit, an individual room works here in the so-called emergency mode, by which the call release is indicated optically outside the room via the indicator lamp and in the rooms with set presence the call is displayed in text form with the call location. System faults are indicated immediately to the care staff as well as to the service personnel via the various display units. Thereby, the system can also be supported through remote maintenance by external service technicians.

#### Installation

The entire installation of the individual system components occurs via a new installation concept by "Plug-and-Play". This makes the whole installation process noticeably easier.

The possibility of mistakes is almost impossible and the simplified maintenance and administration of the system contribute to more efficiency.

In order to keep the installation of the system as simple and economical as possible, the components have been conceived in such a way that all units in the rooms can be installed in connection with round switch boxes. Depending on the various structural characteristics of the building, surface, flush and cavity-wall mounting can be used. The wires are connected by means of screw plug-in terminals.

A particular feature is the construction of the pull cord switch, which can be installed on the wall or the ceiling. This means that it is possible to install the pull cord switch over the lid supports of the WC in the disabled WC.



#### Operation

The call system also proves its worth in operation through its economic advantages. On the one hand, the energy-saving technology of the LED corridor indicator lamps leads to less energy consumption than for conventional lamps, and on the other hand, maintenance expenditure is reduced since the operating life of LED's is up to 20 times longer than that of the previously used piston lamps or soffits.



# **P7 System Overview (with Speech Function)**



# **P7** System Overview (without Speech Function)





# **P7 System Control Unit**

The P7 System Control Unit is the master control unit for the entire system and performs the necessary monitoring functions as well as inter-ward communication. The system control unit establishes the connection to the cross-linked computer architecture via a LAN connection.

This connection enables system configuration and data exchange with optionally available workstations on the wards. If required, data can also be exchanged with existing information systems. On request, appropriate hardware and software interfaces are available for forwarding call information to messaging systems (pagers, wireless DECT telephone systems etc.). System configuration is enabled via LAN connection using standardised browser technology.

The P7 System Control Unit is equipped with modern processor technology and consists of a computer and a fuse holder for top hat rail mounting. No movable storage media are used. The system and operating data are stored in non-volatile memory media in flash technology. Likewise, the system control unit comprises the corresponding software for call and presence logging. The data is stored in the system control unit. Protocol evaluation can be made at any connected duty room monitor or on-site PC via a standard LAN connection. Remote maintenance is possible for the on-site PC via standard LAN connection.

#### The system control unit includes:

- $\hfill\square$  Software for call and presence logging, as well as protocol evaluation
- □ Configuration software for programming the system parameters of the nurse call light system

#### **TECHNICAL DATA:**

| Туре:                   | SSZ-P7-Typ5                                 |
|-------------------------|---|
| Order number:           | 30 01 02                                    |
| Installation:           | Mounted on DIN top hat rail                 |
| Power supply:           | Power max. 10 Watt<br>Input voltage 24 Volt |
| Dimensions (W x H x D): | 45 x 110 x 90 mm                            |

# Mobile call display

Messages can be forwarded to an Android mobile phone for viewing and processing. The messages are also indicated even though the background application or lock screen is enabled. A prerequisite for this function is blanket WLAN coverage throughout the entire facility.

| Order number: | 30 09 06         |
|---------------|------------------|
| Туре:         | SOFT-ANZEIGE-APP |

# **SIP-Gateway**

The SIP gateway connects the call system to the telecommunications system. Calls can be queried and addressed with mobile telephone devices in other associated rooms or by other beds.

#### FEATURES:

| $\square$ SIP gateway for connection to the telecomm | unications system |
|--|-------------------|
| □ Speech channel per ward controller                 |                   |
| Order number:  | 30 09 02          |
| Туре:  | SOFT-SIP-P7       |



# **FAP Interface**

For the system control unit to connect the call system to the fire alarm system using ESPA 4.4.4. protocol.

| Order number: | 30 09 01    |
|---------------|-------------|
| Туре:         | SOFT-BMA-P7 |

# **DECT Interface**

For the system control unit for connection of the call system to the control unit of a DECT system or to the message server of a telecommunications system with integrated DECT devices using ESPA 4.4.4 protocol.

| Order number: | 30 09 00     |
|---------------|--------------|
| Туре:         | SOFT-DECT-P7 |

# Graphic message display

Call and presence display in a ward graphic.

| Order number: | 30 09 07          |
|---------------|-------------------|
| TYPE:         | SOFT-STATIONSPLAN |

# Interface DCF77 ISO

The DCF77 module consists of software and hardware components for coupling the DCF77 signal to the patient call system P7 via an interface on the system central control unit P7.

| Order number: | 39 10 03      |
|---------------|---------------|
| TYPE:         | DCF-SW-HW-ISO |

# **19" Bracket**

With DIN rail for installation of the control technology in 19" cabinet (or similar).

| Order number: | 32 09 00       |
|---------------|----------------|
| Туре:         | 19"-TRÄGER-3HE |

# Isolator

The isolator is used when mounting an additional top hat rail onto the 19-inch bracket (Order number 320900) and intended for electrical isolation of central nurse call light components to DIN 0834 in a 19-inch cabinet.

| Order number: | 32 09 10 |
|---------------|----------|
| Туре:         | ISOLATOR |







# **SBUS Splitter**

The P7-SBus splitter is the solution for star-shaped cabling already in place at the customer's facility. The ward bus can now be divided into up to 8 lines with a maximum of 60 sharing units each.

### FUNCTION:

Cabling of audio & data in the ward data bus, each via separate splitters

For use as a repeater in systems without speech function in case of long distances

Two stage cascading to increase the number of spurs (max. 2 splitters between ward bus unit and ward controller)

#### **TECHNICAL DATA:**

| Connections:            | 10 double-pole screw terminals for<br>power supply, primary connection to<br>the central unit and connection of the<br>spurs |
|-------------------------|--|
| Number of bus sharing u | nits: max. 60  |
| Number of harnesses pe  | r structure: max. 8  |
| Dimensions:             | 121 x 55 x 24 mm   |
| Colour:                 | silver, similar to RAL 9006  |
|                         |  |
| Order number:           | 30 00 30   |
| Туре:                   | SBUS-SPLITTER-P7   |





# RS232 Isolator 5kV

The RS232-ISO converter has been designed to protect against incalculable currents and voltages via the RS232 line of a medical facility. It enables the electrical isolation between devices for medical equipment and systems according to DIN EN 60601-1. The converter can be easily inserted into an existing cabling concept and behaves completely transparently. An external power source is not required for operation.

#### **TECHNICAL DATA:**

| Standards:                           | EIA/TIA RS-232 CCITTV.24 Standard                           |
|--------------------------------------|---|
| Data throughput:                     | 300bps - 912000bps  |
| Power input:                         | < 10mA (drawn at the interface)                             |
| Connections:                         | RS-232 DCE/ RS-232 DTE                                      |
| Electric strength:                   | 5.000Vrms (AC) for longer than 1 min. (increased isolation) |
| ESD resistance:                      | 15kV  |
| Surge protection:                    | 600W/sec  |
| Conforms to the following standards: | IEC / DIN EN 60601-1 /<br>DIN EN 60950-1 + IEC              |

#### **MECHANICAL DETAILS:**

| Dimensions: | 62.8 x 33.8 x 17.8 mm |
|-------------|-----------------------|
| Material:   | PC-ABS                |
| Colour:     | black                 |
|             |                       |

| Order number: | 37 41 21             |  |
|---------------|----------------------|--|
| Туре:         | RS232-ISO-KONVER-5KV |  |







For the power supply of the P3/P7 call systems.

- Closed design
- □ Contact-protected screw connections
- □ Current limiting overload protection, auto recovery, short-circuit proof
- □ Integral noise filter, minimal ripple
- □ Integral earth fault monitoring
- □ Including miniature-fuse holder for top hat rail mounting for secondary protection
- □ Fault / power failure indication in accordance with DIN VDE 0834 Parts 1+2
- □ DIN EN 60601-1/2xMOPP

#### **TECHNICAL DATA:**

| Туре:                   | NG-5A-TYP4                |
|-------------------------|---------------------------|
| Order number:           | 32 04 01                  |
| Weight:                 | 0.92 kg                   |
| Installation:           | Mounted on 35 mm DIN rail |
| Input current:          | 0.7 A (230V AC)           |
| Input voltage:          | 90-265 V AC, 47-63 Hz     |
| Output current:         | 5 A                       |
| Output voltage:         | 24.5 V DC                 |
| Dimensions (W x H x D): | 146 x 205 x 100 mm        |

# Mains Adapter 24V 10A

For the power supply of the P3/P7 call systems

- □ Closed design
- □ Contact-protected screw connections
- □ Current limiting overload protection, auto recovery, short-circuit proof
- □ Integral noise filter, minimal ripple
- □ Integral earth fault monitoring
- □ Including miniature fuse-holder for top hat rail mounting for secondary protection
- □ Fault / power failure indication in accordance with DIN VDE 0834 Parts 1+2
- DIN EN 60601-1/2xMOPP

#### **TECHNICAL DATA:**

| Туре:                   | NG-10A-TYP4               |
|-------------------------|---------------------------|
| Order number:           | 32 04 11                  |
| Weight:                 | 1.56 kg                   |
| Installation:           | Mounted on 35 mm DIN rail |
| Input current:          | 1.3 A (230V AC)           |
| Input voltage:          | 90-265 V AC, 47-63 Hz     |
| Output current:         | 10 A                      |
| Output voltage:         | 24.5 V DC                 |
| Dimensions (W x H x D): | 146 x 205 x 100 mm        |

# **Power Rectifier with UPS Module** 24V, 5A, incl. Rechargeable Battery

For the power supply of the P3/P7 call systems

- □ Rectifier with integrated UPS module comprising:
- □ mains adapter, UPS module
- $\hfill\square$  Ensures uninterrupted output coupling to battery in case of total power failure
- □ With automatic total battery discharge protection
- □ Automatic battery testing
- □ Temperature-controlled battery-charging
- Earth fault monitoring
- □ Including miniature-fuse holder for top hat rail mounting for secondary protection
- □ Fault / power failure indication in accordance with DIN VDE 0834 Parts 1+2
- □ DIN EN 60601-1/2xMOPP

#### **TECHNICAL DATA:**

| Туре:                   | NG-USV-5A-SET4            |
|-------------------------|---------------------------|
| Order number:           | 32 04 49                  |
| Weight:                 | 8.5 kg                    |
| Installation:           | Mounted on 35 mm DIN rail |
| Input current:          | 0.7 A (230V AC)           |
| Input voltage:          | 90-265 V AC, 47-63 Hz     |
| Output current:         | 5 A                       |
| Output voltage:         | 24.5 V DC                 |
| Dimensions (W x H x D): | 146 x 205 x 100 mm        |

# **Rechargeable Battery** for BM-5A-TYP4

| Order number: | 32 04 55        |  |
|---------------|-----------------|--|
| Туре:         | ERSATZAKKU-TYP4 |  |



# **Power Rectifier with UPS Module** 24V, 10A, incl. Rechargeable Battery

For the power supply of the P3/P7 call systems

- □ Rectifier with integrated UPS module comprising:
- □ mains adapter, UPS module
- □ Ensures uninterrupted output coupling to battery in case of total power failure
- □ With automatic total battery discharge protection
- □ Automatic battery testing
- □ Temperature-controlled battery-charging
- □ Earth fault monitoring
- □ Including miniature-fuse holder for top hat rail mounting for secondary protection
- □ Fault / power failure indication in accordance with DIN VDE 0834 Parts 1+2
- □ DIN EN 60601-1/2xMOPP

#### **TECHNICAL DATA:**

| Dimensions (W x H x D): | 146 x 205 x 100 mm        |
|-------------------------|---------------------------|
| Output voltage:         | 24.5 V DC                 |
| Output current:         | 10 A                      |
| Input voltage:          | 90-265 V AC, 47-63 Hz     |
| Input current:          | 1.4 A (230V AC)           |
| Installation:           | Mounted on 35 mm DIN rail |
| Weight:                 | 12 kg                     |
| Order number:           | 32 04 59                  |
| Туре:                   | NG-USV-10A-SET4           |

# **Rechargeable Battery** for BM-10A-TYP4

**Order number:** Type:

## 32 04 56 **ERSATZAKKU-TYP4**

# Ward Controller P7

The ward controller is the central control unit for any ward and regulates optical and acoustic call signalling in the respective rooms. by communicating with the terminals, modules, displays and interfaces via the ward data bus. It checks the function of the connected units and reports faults in devices or lines.

The ward controller supports a ward's care organisation and provides options of dividing it into any number of care groups. Rooms from other wards can also be integrated. The ward controllers are connected to the system control unit via the system bus.

#### Functions:

- Optical and acoustic call processing
- $\hfill\square$  Monitoring connected devices and lines

#### **MECHANICAL DETAILS:**

| Dimensions (W x H x D): | 40 x 110 x 90 mm   |  |
|-------------------------|--|--|
| Power supply:           | Electrical power max. 15 Watt<br>Input voltage 24 V DC   |  |
| Interfaces:             | <ul> <li>Ward bus to the room terminals</li> <li>System bus to central control unit</li> </ul> |  |
| Installation:           | Mounted on DIN top hat rail  |  |
| Order number:           | 30 00 12   |  |
| Туре:                   | SZ-P7-TYP5   |  |

# **Standalone Ward Controller P7**

As described above, but with additional functionality of the system control unit P7. Connecting further ward controllers is excluded.

#### **MECHANICAL DETAILS:**

| Order number: | 30 01 12      |  |
|---------------|---------------|--|
| Туре:         | SZ-P7-AU TYP5 |  |

# **Serial Interface**

For connection of serial interfaces like BMA and DECT, e.g. via ESPA 4.4.4 protocol.

| Order number: | 37 41 22          |
|---------------|-------------------|
| Туре:         | <b>IF-SERIELL</b> |



# **Galvanic Ethernet Isolator**

This device has been designed to protect patients and nursing staff against unpredictable electrical currents via the network wiring of a medical facility. This is a prerequisite for the safe operation of a device in the patient environment. The isolator was developed for galvanic network isolation for medical electrical devices and systems and to meet the requirements of IEC60601-1 and IEC60950-1. The isolation ensures non-restrictive behaviour in the network. The installation is implemented by simply inserting the isolator into an existing network cabling. An external power source is not required. There is no need to install drivers for operation.

#### **TECHNICAL DATA:**

| Data throughput:                     | 10/100/1000 Mbit/S  |
|--------------------------------------|---|
| Electric strength:                   | 4,000Vrms (AC) for longer than 1 min./ 50/60Hz for longer than 1 min. |
| Isolation:                           | Increased isolation at 250 Vrms. operat                               |
| Conforms to the following standards: | IEC / DIN EN 60601-1 / and IEC / DIN E                                |
| MECHANICAL DETAILS:                  |   |
| Dimensions:                          | 72 x 50 x 22 mm   |

| Туре:         |       | ETHERNET |
|---------------|-------|----------|
| Order number: |       | 30 0     |
| Colour:       | black |          |
|               |       |          |

# **Ethernet Switch 8 Port**

For coupling up to 7 units: ward controllers, duty room monitors, duty room consoles or ward PCs already in place for call notification via LAN interface with the system control unit P7.

#### **MECHANICAL DETAILS:**

| Type: SWITCH-8   |
|--|
| Order number: 30 00 2  |
| Colour: dark gre   |
| □ Operating temperature: -40 to 75°C   |
| $\square$ Mounted on DIN rail (top hat rail) or wall mounted                     |
| $\square$ IP 30 protection with metal housing                                    |
| $\hfill\square$ Redundant current supply: 24 V DC with polarity reversal protect |
| B-port 10 / 100 Mbps Auto MDI/MDI-X RJ45   |
| Dimensions (W x H x D): 140 x 95 x 37 mm   |

# **Ethernet Switch 5 Port**

For coupling up to 4 units: ward controllers, duty room monitors, duty room consoles or ward PCs already in place for call notification via LAN interface with the system control unit P7.

#### **MECHANICAL DETAILS:**

| Туре:  | SWITCH-5-F     |
|--|----------------|
| Order number:  | 30 00 20       |
| Colour:  | dark grey      |
| □ Operating temperature: -40 to 75°C                   |                |
| □ Mounted on DIN rail (top hat rail) or wall mounted   |                |
| $\Box$ IP 30 protection with metal housing             |                |
| □ Redundant current supply: 24 V DC with polarity reve | rsal protectio |
| □ 5-port 10 / 100 Mbps Auto MDI/MDI-X RJ45             |                |
| Dimensions (W x H x D): 140 x 95 x 37 mm               |                |



5656V

ting voltage

EN 60950-1

00 22 I-TRENNER





y D •P7



# **Duty Room Monitor**

The hospicall duty room monitor P7 is a 16:9 touch screen PC for extended display and operation of the P7 call system. This extension enables uncomplicated and fast processing of calls and also includes additional functions. The intuitive menu/system guidance ensures rational processes in everyday nursing routine and noticeably reduces the staff's workload.

It is connected to the Ethernet system bus (LAN 1) via a separate bus junction box.



#### FUNCTIONS:

| ~ |
|---|
| ~ |
| ~ |
| ~ |
| ~ |
| ~ |
| ~ |
| ~ |
|   |

#### **MECHANICAL DETAILS:**

| Туре:                   | DZM7 Monitor   |
|-------------------------|--|
| Order number:           | 30 02 06   |
| Housing:                | PC with touch screen, pedestal and connecting cables |
| Material:               | plastic  |
| Protection class:       | Class II/IP65, IPX1 housing                          |
| Dimensions (L x W x H): | 394 x 257 x 57 mm                                    |
|                         |  |



# **Room Terminal Duty Room with Display** with/without Speech Function

The room terminal for duty room is the central control and operating unit for the duty room. Large membrane-covered push-buttons for presence marking as well as for call/emergency call release are available for the care staff. Targeted important information for the care staff is indicated via a backlit, colour, graphic touch-screen display. Thereby, the touch-screen can be used to call up individual additional functions and to link wards linked together. For optical call indication in the corridor, a corridor indicator lamp or a nameplate with integrated corridor indicator lamp can be connected to the duty room terminal via the room bus. Via the ward bus, the duty room terminal of a ward is connected to the associated ward controller. A loudspeaker and microphone enable voice communication with residents, patients and care staff by using the latest DSP technology and digital voice transmission.



Room terminal duty room with display and speech function



Room terminal duty with RT-AW and

speech function



Room terminal duty without speech function

| green        | green   |
|--------------|---|
| green        | green   |
| yellow       | -   |
| red + red    | red + red   |
| blue         | -   |
| $\checkmark$ | $\checkmark$  |
| -            | -   |
| $\checkmark$ | V   |
|              | green<br>green<br>yellow<br>red + red<br>blue<br>V<br>V |

ZT70-D

# FUNCTIONS

| Presence / cancel push-buttons   | green                 |
|--|-----------------------|
| With remind lamp   | green                 |
| Presence push-button with remind lamp                                    | yellow                |
| Call / emergency call push-button with nurse symbol and reassurance lamp | red + red             |
| Doctor call push-button with doctor symbol and reassurance lamp          | blue                  |
| Call indication via colour graphic touch-screen display                  | $\checkmark$          |
| Acoustic call forwarding through noise generator                         | ~                     |
| Emergency mode by failure of data traffic                                | $\checkmark$          |
| Ward linking via graphic touch-screen display                            | $\checkmark$          |
| Loudspeaker and microphone   | ~                     |
| Acoustic and optical indication of faults in room via display            | <ul> <li>✓</li> </ul> |

#### **MECHANICAL DETAILS:**

| Dimensions (L x W x H):   | 150 x 100 x 30 mm  | 150 x 100 x 30 mm  | 150 x 100 x 30 mm   |
|---|--|--|---|
| Protection class:   | IP40   | IP40   | IP40  |
| Material:   | PC-ABS   | PC-ABS   | PC-ABS  |
| Housing:  | white, similar to RAL 9010   | white, similar to RAL 9010   | white, similar to RAL                                     |
| Call push-button:   | red, similar to RAL 3000   | red, similar to RAL 3000   | red, similar to RAL 3                                     |
| Doctor call push-button:  | blue, similar to RAL 5010  | -  | blue, similar to RAL                                      |
| Presence push-button:   | green, similar to RAL 6018   | green, similar to RAL 6018   | green, similar to RAL                                     |
| Presence push-button:   | yellow, similar to RAL 1023  | -  | yellow, similar to RAL                                    |
| Connection:   | <ul> <li>Ward bus to ward controller</li> <li>Room bus</li> <li>Audio bus</li> </ul> | <ul> <li>Ward bus to ward controller</li> <li>Room bus</li> <li>Audio bus</li> </ul> | <ul><li>□ Ward bus to ward o</li><li>□ Room bus</li></ul> |
| Installation:<br>The room terminal is connected to the supply network by means of screw / plu<br>board. | g-in terminals and then clipped into the connection                                  | V  | V   |
| Connection board for duty room terminal included in delivery  | ✓  | V  | v   |
| Order number:   | 30 05 20   | 30 05 21   | 30 05 25  |
| Туре:   | ZT7-D  | 7T7-D2T  | 7T70-D  |

# **Connection Board for Room Terminal**

#### Included in delivery





# **Table Top Stand for Duty Room Terminal**

ZT7-D2T

The table top stand with a 2 m connection cable is intended to accommodate the ward's duty room terminal which is plugged into the stand and then locked into place. A separate bus junction box serves to connect the duty room terminal with table top stand with the ward bus.



#### Room terminal duty with RT-AW without speech function

| x 30 mm                    | 150 x 100 x 30 mm   |
|----------------------------|---|
|                            | IP40  |
|                            | PC-ABS  |
| ilar to RAL 9010           | white, similar to RAL 9010  |
| r to RAL 3000              | red, similar to RAL 3000  |
| ar to RAL 5010             | -   |
| ilar to RAL 6018           | green, similar to RAL 6018  |
| ilar to RAL 1023           | -   |
| s to ward controller<br>us | <ul> <li>Ward bus to ward controller</li> <li>Room bus</li> </ul> |

| ~       |    |
|---------|----|
| ~       |    |
| 30 05 2 | 26 |
| ZT70-D2 | 21 |

| Гуре:                   | TISCHPULT P7                              |
|-------------------------|---|
| Order number:           | 30 02 15                                  |
| Installation:           | via plug-in connector in bus junction box |
| Colour:                 | milky white, similar to RAL 9010          |
| Protection class:       | Class IP40                                |
| Dimensions (L x W x H): | 200 x 120 x 100 mm                        |
|                         |   |

# **Room Terminal with Display with / without Speech Function**

The room terminal is the central control and operating unit for the patient / resident room. Large membrane-covered push-buttons with clear pressure points for presence marking as well as for call/emergency call release provide a good tactile feel. The colour display of the room terminal is supplied with a touchscreen function, so that speech communication can be set up with the patient or person releasing the call by touching the displayed call. The most important information is written large, so that it is easily readable from a distance and easy and secure operation is guaranteed. For optical indication of a call in the corridor, a corridor indicator lamp or a nameplate with integrated corridor indicator lamp can be connected to the room terminal via the room bus. Via the ward bus, the room terminal of a ward is connected to the associated ward controller. A loudspeaker and microphone enable voice communication with residents, patients and care staff by using the latest DSP technology and digital voice transmission.









Room terminal P7 with RT-AW

**Room terminal P7 without speech** function

#### FUNCTIONS:

| Presence / cancel push-buttons   | green                | green     | green     | green        |
|--|----------------------|-----------|-----------|--------------|
| Remind lamp  | green                | green     | green     | green        |
| Presence push-button with remind lamp                                    | yellow               | -         | yellow    | -            |
| Call / emergency call push-button with nurse symbol and reassurance lamp | red + red            | red + red | red + red | red + red    |
| Doctor call push-button with doctor symbol and reassurance lamp          | blue                 | -         | blue      | -            |
| Call indication via colour graphic touch-screen display                  | <ul> <li></li> </ul> | v         | v         | $\checkmark$ |
| Acoustic call forwarding through noise generator                         | <ul> <li></li> </ul> | v         | v         | $\checkmark$ |
| Emergency mode by failure of data traffic                                | <ul> <li></li> </ul> | v         | v         | $\checkmark$ |
| Control of up to 16 active room bus clients                              | <ul> <li></li> </ul> | v         | V         | $\checkmark$ |
| Monitoring of connected call release elements                            | <ul> <li></li> </ul> | v         | v         | ~            |
| Acoustic and optical indication of faults in room via display            | $\checkmark$         | v         | v         | V            |
| Loudspeaker and microphone   | V                    | V         |           | -            |

#### **MECHANICAL DETAILS:**

| Dimensions (L x W x H):  | 150 x 100 x 30 mm  | 150 x 100 x 30 mm  | 150 x 100 x 30 mm   |
|--|--|--|---|
| Protection class:  | IP40   | IP40   | IP40  |
| Material:  | PC-ABS   | PC-ABS   | PC-ABS  |
| Housing:   | white, similar to RAL 9010   | white, similar to RAL 9010   | white, similar to RAL   |
| Call push-button:  | red, similar to RAL 3000   | red, similar to RAL 3000   | red, similar to RAL 3   |
| Doctor call push-button:   | blue, similar to RAL 5010  | -  | blue, similar to RAL S  |
| Presence push-button:  | green, similar to RAL 6018   | green, similar to RAL 6018   | green, similar to RAL   |
| Presence push-button:  | yellow, similar to RAL 1023  | -  | yellow, similar to RAL  |
| Connection:  | <ul> <li>Ward bus to ward controller</li> <li>Room bus to the call units<br/>(max. 16) in room</li> <li>Audio bus</li> </ul> | <ul> <li>Ward bus to ward controller</li> <li>Room bus to the call units<br/>(max. 16) in room</li> <li>Audio bus</li> </ul> | <ul> <li>Ward bus to ward of</li> <li>Room bus to the of (max. 16) in room</li> </ul> |
| Installation:<br>The room terminal is connected to the supply network by means of screw / pl<br>board. | ug-in terminals and then clipped into the connection   | V  | v   |
| Connection board for duty room terminal included in delivery.  | V  | <b>v</b>   | ~   |
| Order number:  | 30 05 10   | 30 05 11   | 30 05 15  |
| Туре:  | ZT7  | ZT7-2T   | ZT70  |
|  |  |  |   |



Room terminal P7 without speech function with RT-AW



# **Room Modules**

The room module is the central control unit for a room. The operating elements of a room are connected via the room module. Via the ward bus, the room modules are connected to a ward with the corresponding ward controller, as well as with each other. Mixed operation with room terminals and other room modules is possible.



Room module with 4-way LED indicator lamp and bed recognition



Room module with call and presence push-buttons



Room module with call and presence push-buttons and display

red

red

green

green ~ ~

~



### FUNCTIONS:

### **MECHANICAL DETAILS:**

Dimensions: Protection class Colour: Installation:

Order number Type:

| Display of call / emergency call / WC call                                | red/white    | -            |
|---|--------------|--------------|
| Display of presence   | green/yellow | -            |
| Call / emergency call release via call push-<br>button with nurse symbol  | -            | red          |
| Orientation / reassurance lamp  | -            | red          |
| Call cancellation / presence marking via cancel / presence push-button    | -            | green        |
| Remind / presence lamp  | -            | green        |
| Emergency mode by failure of data traffic                                 | $\checkmark$ | V            |
| Acoustic call forwarding through noise generator                          | $\checkmark$ | <b>v</b>     |
| Indication of faults in the room  | $\checkmark$ | $\checkmark$ |
| Monitoring of connected call release elements                             | ~            | <i>v</i>     |
| Display 2 x 8 characters for indication of call<br>location and call type | -            | -            |

#### **MECHANICAL DETAILS:**

| Туре:                                       | ZM71                                   | ZM72   | ZM72-DP                                |
|---|--|--|--|
| Order number:                               | 30 05 30                               | 30 05 40                                     | 30 05 45                               |
| Delivery including cover plate and frame.   | $\checkmark$                           | $\checkmark$                                 | <i>v</i>                               |
| Installation:                               | Mounted in fm switch box DIN 4907      | 3. All connection terminals are colour-coded | plug terminals.                        |
| Inputs:                                     | Active call units (max 16) on room bus | Active call units (max 16) on room bus       | Active call units (max 16) on room bus |
| Outputs:                                    | LED lamps                              | Buzzer, reassurance lamps, room bus          | Buzzer, reassurance lamps, room bus    |
| Interfaces:                                 | Ward bus for ward controller           | Ward bus for ward controller                 | Ward bus for ward controller           |
| Doctor call push-button with doctor symbol: | -                                      | -  | -                                      |
| Presence 2:                                 | -                                      | -  | -                                      |
| Cancel / presence push-button:              | -                                      | green, similar to RAL 6018                   | green, similar to RAL 6018             |
| Call push-button with nurse symbol:         | -                                      | red, similar to RAL 3000                     | red, similar to RAL 3000               |
| Colour:                                     | white, similar to RAL 9010             | white, similar to RAL 9010                   | white, similar to RAL 9010             |
| Material:                                   | PC-ABS                                 | PC-ABS                                       | PC-ABS                                 |
| Protection class:                           | IP40                                   | IP40   | IP40                                   |
| Dimensions (L x W x H):                     | 81 x 81 x 15 mm                        | 81 x 81 x 15 mm                              | 81 x 81 x 15 mm                        |

# **Room Module** with Blank Cover

The operating elements of a room are connected via the room module. Via the ward bus, a ward's individual room modules are connected both with the corresponding ward controller and with each other.

Mixed operation with room terminals and other room modules is possible.

- Emergency mode in case of data traffic failure
- Monitoring connected call activation elements
- Indication of faults in the room
- Interfaces: Ward bus to ward controller
- Inputs: Active call units (max. 16) in the room bus
- All connection terminals are of the colour-coded plug type

|    | ZM71-BLIND  |
|----|---|
| :  | 30 05 30-Blind  |
|    | Mounted in fm switch box DIN 49073. All connection terminals are colour-coded plug terminals. |
|    | white, similar to RAL 9010  |
| s: | IP40  |
|    | 81 x 81 x 15 mm   |
|    |   |

# **Call Units**

FUNCTIONS:

With call push-buttons in membrane keypad.

Call / emergency call release via call pushbutton with nurse symbol  $% \left( {{\left[ {{{\rm{c}}} \right]}_{{\rm{c}}}}} \right)$ 

Call cancellation / presence marking via cancel

Acoustic call forwarding through noise generator

Socket for call  $\slash$  emergency call release via the attached operating device

Orientation / reassurance lamp

Connection occurs via room bus

/ presence push-button Remind / presence lamp



red

red

~

Call unit



red

red

green

green

~

Call and cancel unit





Call and cancel unit for separate WC call cancellation (acknowledgement only)

Call unit with auxiliary socket

# red red red red green --~ ~ ~

#### **MECHANICAL DETAILS:**

| Dimensions (L x W x H):                   | 81 x 81 x 15 mm            | 81 x 81 x 15 mm            |
|---|----------------------------|----------------------------|
| Protection class:                         | IP40                       | IP40                       |
| Material:                                 | PC-ABS                     | PC-ABS                     |
| Colour (cover plate and frame):           | white, similar to RAL 9010 | white, similar to RAL 9010 |
| Call push-button with nurse symbol:       | red, similar to RAL 3000   | red, similar to RAL 3000   |
| Cancel push-button:                       | -                          | green, similar to RAL 6018 |
| Interfaces:                               | -                          | -                          |
| Installation: in fm switch box, DIN 49073 | ~                          | <b>v</b>                   |
| Delivery including cover plate and frame  | V                          | ~                          |
| Order number:                             | 30 06 10                   | 30 06 00                   |
| Туре:                                     | RT-P7                      | RT-AW-P7                   |

| <u> </u>                   | o. o   | 04 04 4F   |
|----------------------------|--|--|
| 81 X 81 X 15 MM            | 81 X 81 X 15 mm  | 81 X 81 X 15 mm  |
| IP40                       | IP40   | IP40   |
| PC-ABS                     | PC-ABS   | PC-ABS   |
| white, similar to RAL 9010 | white, similar to RAL 9010   | white, similar to RAL 9010   |
| red, similar to RAL 3000   | red, similar to RAL 3000   | red, similar to RAL 3000   |
| green, similar to RAL 6018 | -  | green, similar to RAL 6018   |
| -                          | RJ-45 jack socket for connection of pear<br>and multiple push-buttons or other portable<br>operating devices | RJ-45 jack socket for connection of pear<br>and multiple push-buttons or other portable<br>operating devices |
| ~                          | <i>v</i>   | V  |
| ~                          | V  | V  |
| 30 06 40                   | 30 05 70   | 30 05 75   |
| RT-AT-P7                   | RT-NS-P7   | RT-AT-NS-P7  |



# **Display Module**

The P7 display module is equipped with a red call push-button and a green presence push-button in membrane keypad, as well as a double-spaced display. The module is connected to the room bus by means of colour-coded screw/plug terminals.

#### FUNCTIONS:

Call / emergency call activation via call push-button Orientation / reassurance lamp, red Call cancellation / presence marking via presence push-button Reminder / presence lamp, green Display, 2x8 digits to indicate call type and origin

Acoustic call forwarding via tone generator

Interfaces: room bus to room module ZM71/ZM72







#### Call and cancel unit for separate call cancellation (acknowledgement only)

| red   |  |
|-------|--|
| red   |  |
| green |  |
| -     |  |
| -     |  |
| V     |  |
| V     |  |

| MECHANICAL DET     | AILS:  |
|--------------------|--|
| Dimensions:        | 81 x 81 x 15 mm                                  |
| Protection class:  | IP40   |
| Material:          | PC-ABS   |
| Colour:            | white cover plate and frame, similar to RAL 9010 |
| Cover plate and fr | ame included in the scope of supply              |
| Installation: M    | ounted in fm switch box DIN 49073.               |
| Order number:      | 30 08 20   |
| Туре:              | DM-P7  |

# **Presence / Cancel Units**

With presence / cancel push-button in membrane keypad.



Presence / cancel call push-button, green



Presence / cancel call push-buttons, green and yellow



Presence / cancel call push-button, yellow

yellow

yellow ~



greer

AT-P7

#### FUNCTIONS:

| Presence marking / call cancellation via<br>presence / cancel push-button | green    | green + yellow |
|---|----------|----------------|
| Remind / presence lamp  | green    | green + yellow |
| Acoustic call forwarding through noise generator                          | ~        | V              |
| WC call cancellation via cancel push-button                               | -        | -              |
| Connection occurs via room bus  | <i>v</i> | $\checkmark$   |

| 81 x 81 x 15 mm             | 81 x 81 x 15 mm      |
|-----------------------------|----------------------|
| IP40                        | IP40                 |
| PC-ABS                      | PC-ABS               |
| white, similar to RAL 9010  | white, similar to RA |
| yellow, similar to RAL 1023 | -                    |
| -                           | green, similar to R  |
| ~                           | v                    |
|                             |                      |

#### **MECHANICAL DETAILS:**

| Dimensions (L x W x H):                   | 81 x 81 x 15 mm            | 81 x 81 x 15 mm             |
|---|----------------------------|-----------------------------|
| Protection class:                         | IP40                       | IP40                        |
| Material:                                 | PC-ABS                     | PC-ABS                      |
| Colour (cover plate and frame):           | white, similar to RAL 9010 | white, similar to RAL 9010  |
| Presence / cancel push-button:            | -                          | yellow, similar to RAL 1023 |
| Presence / cancel push-button:            | green, similar to RAL 6018 | green, similar to RAL 6018  |
| Installation: in fm switch box, DIN 49073 | ~                          | ~                           |
| Delivery including cover plate and frame  | $\checkmark$               | V                           |
| Order number:                             | 30 06 50                   | 30 06 70                    |
| Туре:                                     | AW-P7                      | AW-AW2-P7                   |

# **Reader Module**

The reader module has a radio symbol for personification of the nursing staff in patient rooms. Recognition of a transponder is signalled acoustically, i.e. with a single tone for activation and a double tone for cancellation. The device is connected via the room bus.

#### FUNCTIONS:

Call cancellation / presence marking via transponder

Acoustic call forwarding via tone generator

ACCCESSORIES: Order number 322090 Transponder

#### **MECHANICAL DETAILS:**

| Dimensions:       | 81 x 81 x 15 mm                                  |
|-------------------|--|
| Protection class: | IP40   |
| Material:         | PC-ABS   |
| Colour:           | white cover plate and frame, similar to RAL 9010 |



Cover plate and frame included in the scope of supply

| Order number: | 30 07 90 |
|---------------|----------|
| Туре:         | AW-RD-P7 |

# **Transponder for Room Electronics**

30 06 60

AW2--P7

The transponder for "room electronics with reader module" enables contactless activation or cancellation of a set call in a patient room by the nursing staff at a hospital or nursing facility. The transponder works according to the RFID principle.

#### **TECHNICAL DATA:**

| Chip type:          | EM4102/EM4200                         |
|---------------------|---------------------------------------|
| Technology:         | RFID (radio-frequency-identification) |
| Range:              | typically 6 cm                        |
| Frequency:          | 125kHz                                |
| Protection type:    | waterproof IP68                       |
| MECHANICAL DETAILS: |                                       |

| Installation: | Key fob          |
|---------------|------------------|
| Dimensions:   | 45 x 30 x 1.5 mm |
| Colour:       | black            |
| Material:     | Polyamide        |









| Order number: | 32 20 90 |
|---------------|----------|
| Туре:         | TP-ID    |

# **Call Units**

Call units with call push-buttons in membrane keypad.

# **Pull Cord Switches**

# **Pneumatic Call Unit**



Doctor call push-button,

blue



Call push-button red and

doctor call push-button blue



**Push-button combination** 

call push-button, doctor

call push-button,

green

presence yellow and





Pull cord switch is fitted with a red pull cord with knob and a red reassurance / orientation lamp

Pull cord switch with 2 m long pull cord with red knob with nurse symbol and call cancellation via white 80 cm long pull cord with green knob

| -        | -            | -            |
|----------|--------------|--------------|
| red      | red          | red          |
| -        | -            | $\checkmark$ |
| -        | $\checkmark$ | -            |
| -        | -            | -            |
| -        | -            | -            |
| red      | red          | -            |
| -        | -            | red          |
| -        | green        | -            |
| <b>v</b> | $\checkmark$ | $\checkmark$ |
|          |              |              |

## FUNCTIONS:

| Call / emergency call release via call push-button                        | -                | red  | red    |
|---|------------------|------|--------|
| Orientation / reassurance lamp  | -                | red  | red    |
| Doctor call release via doctor call push-button                           | blue             | blue | blue   |
| Presence 1  | -                | -    | green  |
| Presence 2  |                  | -    | yellow |
| Call cancellation / presence marking via<br>cancel / presence push-button | -                | -    | -      |
| Call / emergency call release via pull cord switch                        | -                | -    | -      |
| Call / emergency call release via rubber ball                             | -                | -    | -      |
| Call cancellation via pull cord with knob                                 | -                | -    | -      |
| Connection to room bus via colour-coded screw                             | / plug terminals |      |        |

### **MECHANICAL DETAILS:**

| Dimensions (L x W x H):                     | 81 x 81 x 15 mm            | 81 x 81 x 15 mm            | 81 x 81 x 15 mm             |
|---|----------------------------|----------------------------|-----------------------------|
| Protection class:                           | IP40                       | IP40                       | IP40                        |
| Material:                                   | PC-ABS                     | PC-ABS                     | PC-ABS                      |
| Colour (cover plate and frame):             | white, similar to RAL 9010 | white, similar to RAL 9010 | white, similar to RAL 9010  |
| Call push-button with nurse symbol:         | -                          | red, similar to RAL 3000   | red, similar to RAL 3000    |
| Doctor call push-button with doctor symbol: | blue, similar to RAL 5010  | blue, similar to RAL 5010  | blue, similar to RAL 5010   |
| Presence 1:                                 | -                          | -                          | green, similar to RAL 6018  |
| Presence 2:                                 | -                          | -                          | yellow, similar to RAL 1023 |
| Installation: in fm switch box, DIN 49073   | $\checkmark$               | ~                          | <b>v</b>                    |
| Delivery including cover plate and frame:   | ~                          | V                          | V                           |
| Order number:                               | 30 06 20                   | 30 06 30                   | 30 07 00                    |
| Туре:                                       | DT-P7                      | RT-DT-P7                   | RT-DT-AW-AW2-P7             |

| ZT-P7                       | ZT-AT-P7                    | ZT-P7-PNEU-7                |
|-----------------------------|-----------------------------|-----------------------------|
| 30 06 80                    | 30 06 85                    | 30 06 90                    |
| ✓ incl. flexible cord 2m    | ✓ incl. flexible cord 2m    | V                           |
| ✓ wall and ceiling mounting | ✓ wall and ceiling mounting | ✓ wall and ceiling mounting |
| -                           | -                           | -                           |
| -                           | -                           | -                           |
| -                           | -                           | -                           |
| -                           | -                           | -                           |
| white, similar to RAL 9010  | white, similar to RAL 9010  | white, similar to RAL 9010  |
| PC-ABS                      | PC-ABS                      | PC-ABS                      |
| IP54                        | IP54                        | IP54                        |
| 81 x 81 x 15 mm             | 81 x 81 x 15 mm             | 81 x 81 x 15 mm             |

# **Replacement Pull Cord**

| With red knob with nurse symbol and clip |          |
|--|----------|
| Order number:                            | 31 00 95 |
| Туре:                                    | ZT-SCHNU |

# **Replacement Pull Cord, White**

| With green knob and clip |              |
|--------------------------|--------------|
| Order number:            | 31 00 98     |
| Туре:                    | ZT-SCHNUR-WE |



Pneumatic call unit with pneumatic switch element. Call release through squeezing the red rubber ball on 2 m long clear plastic tube





# **Call Units**

FUNCTIONS:

Call units with call push-buttons in membrane keypad

# **Diagnostic Unit**

Unit with 2 diagnostic sockets 4-pin

Call / emergency call release via call

Socket for call / emergency call release

push-button with nurse symbol Orientation / reassurance lamp

via the attached operating device





~~



Call push-button with auxiliary socket and DIA socket

red

red

~

1

**RT-NS-DIA-P7** 



Call push-button with auxiliary socket and loudspeaker as double unit

81 x 152 x 15 mm

white, similar to RAL 9010

or other portable operating devices

red, similar to RAL 3000

IP40

PC-ABS

| red   | red                   |
|---|-----------------------|
| red   | red                   |
| $\checkmark$  | $\checkmark$          |
| -   | <b>v</b>              |
| $\checkmark$  | <ul> <li>✓</li> </ul> |
| Connection, via colour added carow / plug terminals to the ream bus |                       |

Connection: via colour-coded screw / plug terminals to the room bus

RJ45 socket for connection of pear and multiple push-buttons

~

~

30 04 70

**RT-NS-LS-P7** 

| Connection: via colour-coded screw / plug terminals to the room bus |
|---|
| MECHANICAL DETAILS:   |

Diagnostic socket

Loudspeaker unit

| Dimensions (L x W x H):  | 81 x 81 x 15 mm   | 81 x 81 x 15 mm  |  |
|--|---|--|--|
| Protection class:  | IP40  | IP40   |  |
| Material:  | PC-ABS  | PC-ABS   |  |
| Colour (cover plate and frame):  | white, similar to RAL 9010                                  | white, similar to RAL 9010   |  |
| Call push-button with nurse symbol:  | -   | red, similar to RAL 3000   |  |
| Interfaces:  | 2 x RJ10 sockets for connection of technical medical device | RJ45 socket for connection of pear<br>and multiple push-buttons or other<br>portable operating devices |  |
|  |   | RJ10 socket for connection of technical medical device   |  |
| Installation: in fm switch box, DIN 49073. All connection terminals are colour-coded screw / plug terminalsführt |   |  |  |
| Delivery including cover plate and frame:  | V   | V  |  |
| Order number:  | 30 07 60  | 30 05 80   |  |

| Order number: | 30 07 60 |
|---------------|----------|
| Гуре:         | DIA2-P7  |
|               |          |

# **Loudspeaker Unit**

The loudspeaker unit is used for the free speech communication when the patient handset is in the wall bracket.

The connection is via the audio unit next to it by the bed, e.g. call unit with auxiliary socket.



#### **MECHANICAL DETAILS:**

| Dimensions (L x W x H):   | 81 x 81 x 15 mm            |  |
|---|----------------------------|--|
| Protection class:   | IP40                       |  |
| Material:   | PC-ABS                     |  |
| Colour:   | white, similar to RAL 9010 |  |
| Installation: in fm switch box, DIN 49073 via colour-coded screw / plug terminals |                            |  |
| Order number:   | 30 08 50                   |  |
| Туре:   | LS-P7                      |  |



# **Audio Isolator**

The audio isolator serves to feed the audio signal into the call system. A transformer is used which ensures a reliable potential separation between the audio source / call system in the patient environment.

#### **TECHNICAL DATA:**

| Audio input:           |
|------------------------|
| Conforms to the follow |

#### **MECHANICAL DETAILS:**

| Dimensions:   |  |  |
|---------------|--|--|
| Material:     |  |  |
| Colour:       |  |  |
| Installation: |  |  |
| Order number: |  |  |
| Туре:         |  |  |



Call push-button with auxiliary socket, DIA socket and loudspeaker as double unit



3.5 mm jack bush ing standards: IEC/DIN EN 60601 and IEC/DIN EN 60950-1

> 81 x 81 x 50 mm PC-ABS Cover plate and frame: white, similar to RAL 9010 Mounted in fm switch box DIN 49073. 30 08 55 **AUDIO-ISOLATOR-P7**

# Multiple Push-Button – Tone incl. Self-releasing Protective Adapter

The Multiple Push-Button – Tone is equipped with a call push-button, two light push- buttons, two volume control buttons and a jack socket to connect a head-phone. Included in the supply is a 2 m connection cord and a self-releasing protective adapter with an 8-pole plug for non-destructive plug disconnection.

#### FUNCTIONS:

| Call / emergency call activation via red call push-button with nurse syn | nbo |
|--|-----|
|--|-----|

- Orientation / reassurance lamp, red
- 2 Light push-buttons, yellow
- 2 Push-buttons for volume control, blue

1 Jack socket for headphone connection 3.5 mm

#### **MECHANICAL DETAILS:**

| Туре:                   | MT-TON                     |
|-------------------------|----------------------------|
| Order number:           | 31 04 70                   |
|                         |                            |
| Housing:                | white, similar to RAL 9010 |
| Material:               | ABS                        |
| Protection class:       | IP40                       |
| Dimensions (L x W x H): | 120 x 46 x 22 mm           |



Wall Bracket for

**MECHANICAL DETAILS:** 

Housing:

Type:

**Order number:** 

Multiple Push-Button – Tone

Installation: Screwed or glued connection

Dimensions (L x W x H): 50 x 32 x 23 mm

Wall bracket for Multiple Push-Button & Patient Handset

31 04 47 WH-MT

Also available with service call push-button

white, similar to RAL 9010





#### FUNCTIONS:

Call / emergency call activation via red call push-button

Orientation / reassurance lamp, red

Control of 2 light sources via connected multiple push unit

Audio signal volume control via control push-button of the connected MT-TON

Diagnostic call by activating a potential-free contact in a connected medical device

# Call Unit with Auxiliary Plug Contact and Tone

Functions identical with 300841, but without diagnostic socket

| Order number: | 30 08 40        |  |
|---------------|-----------------|--|
| Type:         | RT-NS-TON-UP-P7 |  |

# Call / Cancel Unit with Auxiliary Plug Contact and Diagnostic Socket

The call / cancel push-button with auxiliary plug contact and diagnostic socket has a red call button in a membrane keypad including a red reassurance lamp, and a green cancel button in a membrane keypad, as well as an 8-pole socket for the connection of pear push-buttons and multiple push-buttons, as well as a 4-pole diagnostic socket.

#### FUNCTIONS:

Call / emergency call activation via red call push-button Orientation / reassurance lamp, red Call cancellation via green cancel push-button Socket for call / emergency call activation via connected consoles Diagnostic call by activating a potential-free contact in a connected medical device

#### INTERFACES:

RJ45 Socket to connect pear and multiple push-button units or other mobile control panels

RJ10 Socket to connect a medical device

#### **MECHANICAL DETAILS:**



Cover plate and frame included in the scope of supply Installation: Mounted in fm switch box DIN 49073

| Order number: | 30 05 85        |
|---------------|-----------------|
| Туре:         | RT-AT-NS-DIA-P7 |

# Call Unit with Auxiliary Plug Contact, Diagnostic Socket and Tone for Trilux Bedside Lights

The call unit with auxiliary plug contact, diagnostic socket and tone is a nurse call light system intended for installation in the Trilux ES 403 lights or Sanesca. It is equipped with a red call button, an 8-pin socket (RJ45) to connect the MT-TON multiple push-button, and a 4-pin socket (RJ10) to connect a medical device, enabling the transmission of an audio signal (TV tone).

#### FUNCTIONS:

| Call / emergency call | activation via red | call push-button |
|-----------------------|--------------------|------------------|
|-----------------------|--------------------|------------------|

Orientation / reassurance lamp, red

Control of 2 light sources via connected multiple push-button

Audio signal volume control via control push-button of the connected  $\ensuremath{\mathsf{MT}}\xspace{\mathsf{TON}}$ 

Diagnostic call by activating a potential-free contact in a connected medical device

34

# Call Unit with Auxiliary Plug Contact, Diagnostic Plug Contact and Tone

The call unit with auxiliary plug contact and tone is equipped with a red call button and an 8-pin socket (RJ45) to connect the multiple push-button MT-TON, as well as a 4-pin diagnostic socket (RJ10), enabling the transmission of an audio signal (TV tone).

#### **MECHANICAL DETAILS:**

| Dimensions (L x W x H): 81 x 81 x 15 mm               |                       |                            |  |
|---|-----------------------|----------------------------|--|
| Protection class:                                     | IP40                  |                            |  |
| Material:   | PC-ABS                |                            |  |
| Colour: v   | white cover plate and | frame, similar to RAL 9010 |  |
| Call push-button:                                     | red with nurse s      | ymbol, similar to RAL 3000 |  |
| Cancel push-butt                                      | on: green, similar to | RAL 6018                   |  |
| Cover plate and frame included in the scope of supply |                       |                            |  |
| Installation: Mou                                     | nted in fm switch box | DIN 49073                  |  |
| Order number:   |                       | 30 08 41                   |  |
| Туре:   |                       | RT-NS-DIA-TON-UP-P7        |  |

# Call Unit with Auxiliary Plug Contact and Tone for Trilux Bedside Lights

Functions identical with 300844, but without diagnostic socket

| Order number: | 30 08 43         |  |
|---------------|------------------|--|
| Туре:         | RT-NS-TON-MVE-P7 |  |



#### **MECHANICAL DETAILS:**

| Туре:                   | RT-NS-DIA-TON-MVE-P7       |
|-------------------------|----------------------------|
| Order number:           | 30 08 44                   |
| Cover plate colour:     | white, similar to RAL 9010 |
| Material:               | PA/PC                      |
| Protection class:       | IP40                       |
| Dimensions (L x W x H): | 80 x 63 x 30 mm            |

# **Patient Handset P7**

Patient handset with auto-release plug system for discreet two-way communication between the patient and the care staff.

# Wall Bracket

In the wall bracket the patient handset is in the parked position, which results in voice communication being switched automatically from discreet to free two-way communication.



# **Pear Push-Button / Multiple Push-Button**



Pear push-button



Multiple push-button with one light button

#### FUNCTIONS:

| Call / emergency call release via call push-<br>button with nurse symbol: | red    | - | red | red    | red    |
|---|--------|---|-----|--------|--------|
| Orientation / reassurance lamp  | red    | - | red | red    | red    |
| Light push-button for switching room light                                | yellow | - | -   | -      | yellow |
| Light push-button for switching reading light                             | yellow | - | -   | yellow | yellow |
| Loudspeaker and microphone for discreet com-<br>munication connection     | ~      | - | -   | -      | -      |

#### **MECHANICAL DETAILS:**

| Dimensions (L x W x H): | 125 x 42 x 20 mm  | 20 x 42 x 29 mm            |
|-------------------------|---|----------------------------|
| Protection class:       | IP40  | IP40                       |
| Material:               | PC-ABS  | PC-ABS                     |
| Cover plate and frame:  | white, similar to RAL 9010  | white, similar to RAL 9010 |
| Installation:           | -   | Screw or glue fixing       |
| Delivery:               | Incl. flexible cord 3m with auto-<br>release protective adapter for<br>non-destructive plug-out |                            |
| Order number:           | 30 05 90  | 30 05 91                   |
| Туре:                   | PHG7  | WH-PHG7                    |

| 110 x 40 x 22 mm           | 110 x 40 x 22 mm           | 110 x 40 x 22 mm           |
|----------------------------|----------------------------|----------------------------|
| IP65                       | IP65                       | IP65                       |
| ABS                        | ABS                        | ABS                        |
| white, similar to RAL 9010 | white, similar to RAL 9010 | white, similar to RAL 9010 |
|                            |                            |                            |

Pear and multiple push-buttons are delivered with 3m long flexible cord and 8-pin auto-release protective adapter for non-destructive plug-out

| 31 04 01 | 31 04 11  |
|----------|-----------|
| BT-hc    | BT-1LT-hc |

# **Connection Extension**

Through the use of connection extensions with 7-pin plug and 7-pin coupling between the connection cord of an operating device and the protective adapter, the connection cords of the operating devices can be extended easily as required.

Call-circuit monitoring enables each disconnection of a connector to be indicated immediately.

#### **MECHANICAL DETAILS:**

| Protection class: | IP 40<br>PV/C              |
|-------------------|----------------------------|
| Colour:           | white, similar to RAL 9010 |
| Order number:     | 31 04 40                   |
| Туре:             | ANSCHL-VERL-2M             |

# Wall Bracket for Pear and Multiple Push-Buttons

| Wall bracket for pear and multiple push-buttons. | Mounting: screw or glue fixing. |
|--|---------------------------------|
| Order number:                                    | 31 04 46                        |
| Туре:  | WH-BT                           |

# **Trapeze Bar Clip**

Trapeze Bar Clip to connect the cable of a patient hand set to the trapeze bar

| Dimensions (L x W x H): | 39 x 66 x 44 mm            |
|-------------------------|----------------------------|
| Colour:                 | white, similar to RAL 9010 |
| Order number:           | 31 04 29                   |
| Туре:                   | BETTGALGEN-C               |
|                         |                            |



#### Multiple push-button with two light buttons

| 31  | 04  | 21   |
|-----|-----|------|
| BT- | 2LT | '-hc |







# **LED Corridor Indicator Lamps**

The corridor indicator lamp has 4 lamp cells, which are fitted with white, red, green and yellow light-emitting diodes. White for indication of bathroom / WC calls.

# **Door Plate with LED Corridor Indicator Lamps**

# Nameplate with LED **Corridor Indicator Lamps**

In a homogeneous design for a homely ambience, with easily changeable nameplates without additional tools and corridor indicator lamps with high-grade LED's.



LED corridor indicator lamp white / red / green / yellow



Door plate with LED corridor indicator lamp white / red / green / yellow



Nameplate with LED corridor indicator lamp white / red / green / yellow

#### FUNCTIONS:

| Display of call / emergency call | red       | red       | red    |
|----------------------------------|-----------|-----------|--------|
| Display of Presence 1            | green     | green     | green  |
| Display of Presence 2            | yellow    | yellow    | yellow |
| Display of bathroom / WC call    | white/red | white/red | white  |

#### **MECHANICAL DETAILS:**

| Dimensions (L x W x H):                   | 81 x 81 x 35 mm                        | 81 x 152 x 35 mm                  |
|---|--|-----------------------------------|
| Protection class:                         | IP40                                   | IP40                              |
| Material:                                 | PC-ABS                                 | PC-ABS                            |
| Colour (cover plate and frame):           | white, similar to RAL 9010             | white, similar to RAL 9010        |
|   |  |                                   |
| Connection:                               |  | To a room terminal or room module |
| Installation:                             | Mounted in fm switch box, DIN<br>49073 |                                   |
| Delivery including cover plate and frame: | v                                      | ~                                 |
| Order number:                             | 30 07 20                               | 30 08 00                          |
| Туре:                                     | ZL-P7                                  | TÜRSCHILD-ZL-P7                   |

| 100 x 148 x 35 mm  |
|--|
| IP40   |
| PC/PA-GF   |
| white, similar to RAL 9010   |
|  |
|  |
| Snap-on fastening on supplied connection board for fixing to fm switch box DIN 49073 |
| -  |
| 30 08 05   |

NAMEN-ZL-4-P7

**Electronics for Connection to a Corridor Indicator Lamp** white / red / green / yellow P7

Adapter for coupling of secondary light sources (24 V DC or 12 V DC) to a room terminal or room module P7. The adapter is supplied with shrinkage tubing and is placed in the fm switch box of the respective secondary light.

#### **MECHANICAL DETAILS:**

| Type:                   | ZL-FL-P7                    |
|-------------------------|-----------------------------|
| Order number:           | 30 07 25                    |
| Installation:           | In fm switch box, DIN 49073 |
| Material:               | PC-ABS                      |
| Dimensions (L x W x H): | 43 x 23 x 22 mm             |
|                         |                             |



# **Nameplate Service Set**

Set consisting of: 100 sheets A4 card, white, fine perforated for easy separation of the nameplate inserts + 10 transparent covers, foiled, software for printing nameplate.

**Order number:** 

Type:



31 07 01 NAMEN-ZL-S

# Corridor Display, Single-Sided for Wall Mounting Corridor Display, Double-Sided for Ceiling Mounting

With dot-matrix LED display for self-illuminating plain text display and integrated noise generator (buzzer).

#### FUNCTIONS DEPENDING ON CONFIGURATION:

Volume level time-related, infinitely variable, configurable Time display, when no calls exist Possibility of colour differentiation of display of calls (red) and time (green)

#### EQUIPPED WITH:



| 13-character LED display (per side)             |
|---|
| Contrast-increased, anti-reflective front panel |

Up to 13-character dot-matrix LED display

Buzzer

With integrated interface for ward bus

#### **MECHANICAL DETAILS:**

| Model:                  | Corridor display, single-sided for wall mounting | Corridor display, double-sided for ceiling mounting                                 |
|-------------------------|--|---|
| Dimensions (L x W x H): | 660 x 100 x 65 mm                                | 660 x 100 x 100 mm  |
| Protection class:       | IP 40  | IP 40   |
| Material:               | steel plate / Al                                 | steel plate / Al  |
| Colour:                 | white, similar to RAL 9010                       | white, similar to RAL 9010  |
| Installation:           | wall mounted                                     | Ceiling mounting, pendulum tubes (length): 500 mm, infinitely adjustable            |
| Accessories:            |  | Including baldachin and fastening material for mounting directly below the ceilings |
| Order number:           | 30 08 70   | 30 08 80  |
| Туре:                   | DP-WAND-P7                                       | DP-DOPPEL-P7  |





# **Interface Door Intercom**

The interface can be used as an independent door intercom or as an interface to a customerinstalled door intercom. Activation occurs via the room bus of an assigned room module.

#### FUNCTION:

- □ Doorbell (only when covered by call push-button)
- □ Microphone amplifier with microphone and loudspeaker
- □ LED as "reassurance lamp"

#### **OUTPUTS:**

- □ Switching of an external LED e.g. as "call recognition"
- $\square$  Relay for switching door opener

#### INPUTS:

- $\hfill\square$  Relay for loudspeaker bus for use as interface
- $\hfill\square$  External microphone input for use as interface

#### **MECHANICAL DETAILS:**

| Dimensions (L x W x D):                  | 81 x 81 x 15 mm            |
|--|----------------------------|
| Protection class:                        | IP40                       |
| Material:                                | PC-ABS                     |
| Colour cover plate and frame:            | white, similar to RAL 9010 |
| Installation:                            | In fm switch box, DIN49073 |
| Delivery including cover plate and frame |                            |
| Order number:                            | 30 07 55                   |
| Туре:                                    | IF-TUER-P7                 |



# I/O-2 Input / Output Module

The interface module has two individually configurable inputs and two individually configurable relay outputs, which can be used in floating or non-floating AC or DC as change-over contacts. Via the opto-decoupled inputs, external alarms (heating, FAP, emergency exit, lift etc.) with a switching voltage of 5 - 30 V AC/DC per contact can be registered in the call system. The relay outputs serve, e.g., to switch external devices (Venetian blinds, lamps, door magnets, horns etc.). The activation of the relay occurs via software according to the configuration in the system control unit. Each output relay can take a load of up to 1A. The switching of the interface module occurs via the room bus.

#### INPUTS:

 $\Box$  2 x opto-decouplers input 5 – 30 V AC/DC

#### OUTPUTS:

 $\Box$  2 x switch contacts max. 1A

 $\square$  2 x LED's for indication of freely programmable applications

#### **MECHANICAL DETAILS:**

| Type                                     |              |
|--|--------------|
| Order number:                            | 30           |
| Delivery including cover plate and frame |              |
| Installation:                            | In fm switch |
| Colour cover plate and frame:            | white, simil |
| Material:                                | PC-ABS       |
| Protection class:                        | IP40         |
| Dimensions (L x W x D):                  | 81 x 81 x 1  |

# I/O-2-AT Input / Output Module with Cancel Push-Button

The module has the same characteristics as the standard version of I/O-2.

The module has an additional green cancel push-button for the cancellation of incoming messages locally.

| Order number: | 30 ( |
|---------------|------|
| Type:         | 102- |

# **Current Impulse Switch**

The current impulse switches are used for switching the room lighting or the reading light. The devices have their own galvanic separated current supply. Installation is either in a double flush mounted switch box or in medical supply units.

#### **TECHNICAL DATA:**

| Dimensions (L x W x H): | 85 x 42 x 32 mm  |
|-------------------------|------------------|
| Operating voltage:      | 150-250 V, 50 Hz |
| Nominal switch rating:  | 6 A / 250 V AC   |

| 04 80  |
|--------|
| [ST-01 |
|        |

| Current impulse switch with two contacts |          |
|--|----------|
| Order number:                            | 32 04 82 |
| Туре:                                    | STST-02  |







L5 mm

lar to RAL 9010 h box, DIN49073

07 50 02-P7

07 51 -AT-P7



# I/O-8 Input / Output Module

The interface module has eight individually configurable inputs and eight individually configurable relay outputs, which can be used in floating or non-floating AC or DC as change-over contacts.

Via the opto-decoupled inputs, external alarms (heating, FAP, emergency exit, lift etc.) with a switching voltage of 5 - 30 V AC/DC per contact can be registered in the call system. The relay outputs serve, e.g., to switch external devices (Venetian blinds, lamps, door magnets, horns etc.). The activation of the relay occurs via software according to the configuration in the system control unit. Each output relay can take a load of up to 8A. The switching of the interface occurs via the ward bus.

#### **MECHANICAL DETAILS:**

| Туре:                         | I08-P7                    |
|-------------------------------|---------------------------|
| Order number:                 | 30 05 60                  |
|                               |                           |
| Installation:                 | Mounted on DIN rail       |
| Colour cover plate and frame: | grey, similar to RAL 7035 |
| Material:                     | PC-ABS                    |
| Current supply:               | via 24 V DC               |
| Protection class:             | IP40                      |
| Dimensions (L x W x H):       | 156 x 90 x 59 mm          |





# **Radio Sets**

The radio sets consist of a pluggable receiver and a mobile transmitter for call release. The transmitter is available as a pendant transmitter or a wrist strap transmitter. The pluggable receiver is connected to a unit with an 8-pin auxiliary socket. A maximum of 32 transmitters can be assigned to a receiver. Frequency: 868.3 MHz

#### Note:

The radio set is an add-on to the call system P7. It does not fulfil the safety requirements of DIN VDE 0834.

# **Radio Set with Wrist Strap Transmitter**

Comprising a wrist strap transmitter and a radio receiver 2 for connection to the 4-pin DIA socket.

| Order number: | 31 20 80    |
|---------------|-------------|
| Туре:         | FUFI2-RJ10- |

# **Radio Set with Pendant Transmitter**

| Order number: | 31 20 10 |
|---------------|----------|
| Туре:         | FUF12    |
|               |          |

# **Radio Set with Wrist Strap Transmitter**

| 31 20 50  |
|-----------|
| FUFI2-AMB |
|           |



# **Radio Receiver 2**

For connection to auxiliary socket P7 8-pin RJ-45

| Order number: | 31 20 20 |
|---------------|----------|
| Туре:         | FU2-STEU |

# **Radio Receiver RJ-10**

For connection to DIA socket P7 4-pin RJ-10

| Order number: | 31 20 60      |
|---------------|---------------|
| Туре:         | FU2-STEU-RJ10 |



# **Radio Pendant Transmitter 2 with Strangulation Protection**

| Order number: | 31 20 30 |
|---------------|----------|
| Туре:         | FU2-UMH  |



| Order number: | 31 20 40 |
|---------------|----------|
| Туре:         | FU2-AMB  |

# **Battery Replacement Set for Radio Transmitter**

Consisting of: special key, sealing ring, 3 V battery (CR2032)

| Order number: | 31 01 95  |
|---------------|-----------|
| Туре:         | FU-FI-BWS |

Warning! The radio transmitters of Series 2 are not compatible with older radio sets (different frequency).







connection to the 4-pin DIA socket.

**Order number:** 

Type:







# **Radio Set with Pendant Transmitter**

Comprising a pendant transmitter with strangulation protection, and a radio receiver 2 for

31 20 70 FUFI2-RJ10-UHS



# **Sound Detector**

The sound detector can initiate a call by directly pressing the push-button as well as via a corresponding noise level. The sensitivity is adjustable. The detector is fitted with a 2 m long connection cord and an RJ-45 connector.

## **EQUIPPED WITH:**

Call / emergency call release via noise level or red call push-button with nurse symbol Orientation / reassurance lamp, red Function indicator lamp, green C

| Order number: | 31 04 60 |
|---------------|----------|
| Туре:         | SW-1     |
|               |          |

# **Sealing Flange for FM Inserts**

With the sealing flange, the increased protection class IP44 can be achieved when using the fm inserts designed for wet rooms. The flange is inserted into the fm box underneath the insert between the support ring and the wall. This allows the moisture from above entering behind the frame to bypass the insert so that it can escape at the bottom without contacting electrical connections or components.

#### FUNCTION:

Type:

| Order number:  | 30 02 60        |
|--|-----------------|
|  |                 |
| Colour:  | transparent     |
| Material:  | PE              |
| Protection class:  | IP44            |
| Dimensions (L x W x H):                                      | 78 x 80 x 25 mm |
| MECHANICAL DETAILS:  |                 |
| Pre-punched openings for screw fastening of the fm socket    |                 |
| Flexible material for easy installation around the fm insert |                 |
| Protection of fm inserts in wet rooms                        |                 |
| For devices with covers larger than 50 x 50 mm               |                 |
| Compliance with the protection class in case of uneven walls |                 |

T

DF-UP

# FUNCTIONS:

2x RJ45 Socket to connect the duty room monitor and console for the room terminal duty room

## MECHAN

Dimensions (L x W x H): 81 x 81 mm Protection class: IP40 Material: PC-ABS Cover plate and frame included in the scope of supply: white, similar to RAL 9010

Order n Type:







# **BUS Junction Socket**

Both the duty room monitor and the console for the room terminal duty room are connected to the ward bus and the system bus via the bus junction socket featuring the 2 x 8 pole connector.

| NICAL | DETAILS: |
|-------|----------|
|       |          |

Mounted in fm switch box, DIN 49073

| number: | 30 02 10 |
|---------|----------|
|         | AD-DZT7  |

# **Design Range**

All push-buttons are supplied complete with cover plates with Other surface designs and colour options ar membrane push-buttons and cover frames in the AS500 ANTI-BAKTERIELL design from Jung, in white, similar to RAL 9010. All room terminals and nameplates with corridor indicator lamps are Please contact us for prices. also typically supplied in white, similar to RAL 9010.

DUNG

DUNG

DUNG

DUNG

DUNG

**GIRA** 

**GIRA** 

**GIRA** 

**GIRA** 

**GIRA** 

A flow

AS500

AS500

**A Creation** 

A550

Standard 55

Standard 55

E2

E2

Event

at an extra charge!

white

white

ivory

white

white

pure white, glossy

pure white, matt

pure white, glossy

pure white, matt

pure white, glossy

| re available on request |                            | B,<br>Berker       | S.1              | polar white, glossy  | XXX XXX-503 |
|-------------------------|----------------------------|--------------------|------------------|----------------------|-------------|
|                         |                            | <b>B</b><br>Berker | S.1              | white glossy         | XXX XXX-504 |
| XXX-XXX-215             |                            | BUSCH-JAEGER       | Busch-balance SI | alpine white         | XXX XXX-105 |
| XXX XXX-201€            |                            | merten             | 1-M              | active white, glossy | XXX XXX-407 |
| XXX XXX-202             |                            | merten             | 1-M              | polar white          | XXX XXX-406 |
| XXX XXX-216             |                            | merten             | M-Smart          | active white, glossy | XXX XXX-405 |
| XXX XXX-217             |                            | merten             | M-Smart          | polar white, glossy  | XXX XXX-401 |
| XXX XXX-301             |                            | merten             | M-Pure           | active white         | XXX XXX-408 |
| XXX XXX-303             |                            | merten             | M-Pure           | polar white          | XXX XXX-409 |
| XXX XXX-305             |                            |                    |                  |                      |             |
| XXX XXX-304             |                            |                    |                  |                      |             |
| XXX XXX-315             | We will gladly assist you. |                    |                  |                      |             |
|                         | P +49 2261 95              |                    |                  |                      |             |



# 

# ... just press the button – and help is on the way!

hospicall GmbH • Martin-Siebert-Straße 1 • 51647 Gummersbach • Germany • T +49 2261 9583 - 100 • info@hospicall.com • www.hospicall.com